

hygea d3ntal.

Ultrawave's **hygea d3ntal** offers outstanding performance for use in dental applications, chiropody, veterinary and many other professional environments, fully compliant to HTM01-05 and ISO standards, to provide reassurance for any professional use.

The **hygea d3ntal** incorporates Ultrawave's *Frequency Leap* Technology, homogenous cavitation activity is generated throughout the ultrasonic bath reducing dead spots and standing waves. This provides a more consistent scrubbing action for excellent results.



The **hygea d3ntal** is a 10 Litre heated, digitally controlled ultrasonic bath.

The menu driven control panel allows simple adjustment of cycle time (0 to 99 minutes), temperature (ambient to 50°C), ultrasonic power and frequency adjustment for precise and accurate processing.

Combining features such as degas mode, fluid level sensor and a memory card port, the **hygea d3ntal** offers repeatable and traceable cleaning performance with consistent cleaning results every time.



Eastgate Business Park
Wentloog Avenue
Cardiff CF3 2EY
United Kingdom

Tel: +44 (0) 845 330 4236
Fax: +44 (0) 845 330 4231
Email: admin@ultrawave.co.uk
Web: www.ultrawave.co.uk

hygea d3ntal

Construction Case	Powder coated 304 Stainless Steel
Construction Tank	304 Stainless Steel
Front Panel	FR ABS Plastic
Control	Digital LCD
Overall Dimensions	345 x 285 x 265mm (WxDxH)
Tank Dimensions	295 x 245 x 150mm (WxDxH)
Basket Dimensions	265 x 215 x 140mm (WxDxH)
Maximum Capacity	10Litres
Working Capacity	9.5Litres
Ultrasonic Power	100W rms (280W p to p)
Ultrasonic Power (per Litre)	11W/Litre
Operating Frequency	32 - 38KHz <i>Frequency Leap</i>
Maximum Heating Capacity	Ambient - 50°C
Heater Power	300W
Timer	0 - 99 minutes
Mains Connection	Standard 3 pin socket, 13A
Drain Outlet	½" BSP VALVE
Lid	ABS Plastic
Basket	Stainless Steel
Gross Weight (inc. lid, basket & packaging)	11.5kg
Net Weight	9kg
Warranty	3 Years

